

Autoclave Report

CE Tank Cleaner™

International Article Number

162234

Product classification

Above the Floor Tool

Composition

Metal Parts: Stainless Steel

Test Conditions

Autoclave

- MLS Series 2.6 cu.ft. Portable Top-Loading Autoclave MLS-3781L-PA

Test Parameters

- 134 °C for 5 min
- Pre-vacuum: 20 kPa
- Cooling-down of product between autoclaving cycles



RESULT

After 10 autoclave cycles the tool does not show any damages or loss of function, but does show some characteristic changes: (See Below)

Hardware:

The hardware experienced minor loosening after autoclave cycles.

Date: 11/30/2018

All information contained herein is given according to our specifications and best knowledge without any warranty or guarantee and can be changed without further notice. The suitability of our products to our customers specific applications and conditions of use has to be determined by our customers. In particular, product users shall not be released from their duty to check all health, safety and environment relevant properties of the delivered goods under their specific conditions of use.